(Keraring

REV 15 - 01/2025

C Keraring

PT-BR	3
PT-PT	8
EN	14
LT	20
ES	25
DE	31
FR	37
IT	43
DA	49
NL	55
EL	61
RO	68
SK	73
SL	79
CZ	84
ET	90
LV	96

SYBOLOGY / SIMBOLOGIA / IMBOLOGIA / SYMBOLIK / SIMBOLOGIA / SYMBOLOGIE / SIMBOLOGIA / SYMBOLOGI / SYMBOLOGIE / ΣΥΜΒΟΛΟΓΙΑ / SIMBOLOGIE / SYMBOLÓGIA / SIMBOLOGIJA / SIMBOLOGIJA / SIMBOLOGIJA / SÛMBOLOGIA / SYMBOLOGIE	101
--	------------

PT **Fabricado por** / EN **Manufactured by** / DE **Hergestellt von** / ES **Fabricado por** /
FR **Fabriqué par** / IT **Fabbricato da** / DA **Fremstillet af** / NL **Gemaakt door** /
EL **Κατασκευασμένο από** / RO **Produs de** / SK **Vyrába** / SL **Proizvajalec** /
LT **Gamintojas** / LV **Ražots** / ET **Toodetud** / CZ **Vyrobeno**:

Mediphacos Indústrias Médicas S/A
Av. Deputado Cristovam Chiarádia, 777. Belo Horizonte, MG
CEP: 30575-815 Brasil – www.mediphacos.com
Tel.: +55 31 2102-2211
CNPJ: 21998885/0001-30
Responsável Técnico: Aline L. Freitas – CRFMG: 23364
ANVISA/MS: 10161000006

MEDIPHACOS

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

WARNING! This device must only be implanted by an experienced eye surgeon trained in the implantation technique of this specific device.

1. DESCRIPTION

Keraring is an Intraströmral Corneal Ring Segment (ICRS, device, implant) are precision implantable devices implanted in the human corneal stroma for the correction of morphological and refractive disorders. The device is manufactured from polymethylmethacrylate (PMMA), and is implanted in the corneal stroma as per specific surgical technique. The device acts upon the corneal tissue by altering its central curvature and shape, thus reducing or eliminating morphological irregularities and astigmatism. The device is composed of one or two semicircular segments of variable arc lengths, apical diameters and thickness. Each ring segment tip has a hole to facilitate manipulation and implantation.

MODEL	DESIGN	SPECIFICATION	DIFFERENTIATION
Keraring SIS e SI6		Intraströmral ring for regularization of corneal ectasia. Arc length: from 60° to 355° Thicknesses: from 0.100 mm to 0.400 mm (from 100 µm to 400 µm) Base width: 0.600mm to 0.800mm	Standard product
Keraring ASS e AS6		Asymmetrical intraströmral ring for regularization of corneal ectasia. Arc: from 90° to 355° Thickness: from 0.100mm to 0.400mm and from 0.400mm to 0.100 mm (from 100 µm to 400 µm and from 400 µm to 100 µm) Base width: 0.600mm to 0.800mm	Variable thickness along the arc. Clockwise (W), Counter-Clockwise (C) and Middle (MT) thickness increase.

MODEL	DESIGN	SPECIFICATION	DIFFERENTIATION
Keraring SG		Intrastromal ring for changing corneal curvature. Arc: from 90° to 340° Curvatures: 22 D to 67 D Thicknesses: from 0.100 mm to 0.400 mm (from 100 µm to 400 µm)	Scheduled change of corneal curvature.

2. INTRODUCTION:

Each Keraring Intrastromal Corneal Ring Segment is individually presented with ethylene oxide sterilized in a neutral plastic case sealed in a sterilization envelope. Implant identification labels are provided and contain model information, dimensions, serial number, lot and sterilization validity. An implant card is also provided.

The product is for single use and should not be re-sterilized or reused. Resterilization using physical (dry heat or vapor) or chemical (ethylene oxide) methods and / or reuse of the product will compromise the physical condition and / or functional characteristics of the device causing harm to the patient.

3. INTENDED PURPOSE

Keraring is an Intrastromal Corneal Ring indicated for Keratoconus treatment in the pediatric (8 years old or older) and adult populations. To be used by ophthalmologists surgeons trained in the Intrastromal Corneal Ring surgical technique. The principle of action of the Intrastromal Corneal Ring is to add mass and volume to the corneal tissue, to reshape its central region towards a flatter and more regular surface. The implantation of the Keraring may be considered in the following cases: Keratoconus in contact lens intolerant patients; Astigmatism following penetrating keratoplasty; Corneal ectasia following Lasik; Pellucid marginal degeneration.

4. INDICATIONS FOR USE

The correct indication for implantation of the Keraring requires a thorough evaluation of topography and pachymetry cornea's conditions, in addition to a complete ophthalmologic examination. The implantation of the device may be considered in the following cases:

- Keratoconus in contact lens intolerant patients;
- Astigmatism following penetrating keratoplasty;
- Corneal ectasia following Lasik;
- Pellucid marginal degeneration

5. CONTRAINDICATIONS

Implantation of Keraring must not be considered in the following conditions:

- Severe keratoconus (keratometry reading steeper than 70.0 D);
- Severe central corneal opacities;
- Hydrops;
- Following penetrating keratoplasty when graft is decentered;
- Severe atopic disease associated with chronic eye rubbing habits;
- Presence of active infection, localized or systemic;
- Autoimmune or immunodeficiency disease;
- Recurrent corneal erosion syndrome;
- Extensive corneal scarring;
- Corneal dystrophy;
- Corneas too thin;
- Vernal Keratoconjunctivitis associated with chronic eye rubbing habits.

NOTE: Mediphacos does not recommend the use of the device for situations that are not indicated in this instruction, as this

may impact product performance and/or generate potential adverse events. In these cases the use will be the responsibility of the user.

6. INTENDED USER

Medical doctors (ophthalmologists) which have knowledge in the ICRS surgical technique.

7. INTENDED PATIENT POPULATION

Keraring can be implanted in children, 8 years old or older, and adults, men and women.

8. PERFORMANCE CHARACTERISTICS

Keraring add mass and volume to the corneal tissue. The device acts upon the corneal tissue by altering its central curvature and shape, thus reducing or eliminating morphological irregularities and astigmatism.

9. EXPECTED CLINICAL BENEFITS

- Treatment of Keratoconus to postpone or avoid invasive procedures;
- Improve visual acuity, quality of vision and quality of life and/or tolerability of glasses or contact lenses;
- Reduces refractive error, cylinder, topographic and keratometer parameters.

10. WARNINGS

- During clinical studies, the following complications were identified: segment extrusion; segment migration; vascularization, opacification and scarring of the tunnel associated with segment' migration; decentration; difficult for segment insertion or for implant; decentration or decentered created tunnel; incomplete tunnel creation / incomplete channel formation; patients unsatisfied because of poor refractive outcome; poor visual outcome; anterior chamber perforation; anterior corneal perforation; corneal perforation; localized infiltrate; misdirection of ring segment in to the stroma; inverted implanted ring; steroid-induced glaucoma; deep implantation; shallow or superficial ring segment placement; asymmetrical positioning of the ring; wrong segment position or segment movement after the procedure; corneal vascularization; corneal neovascularization; night-vision disturbances; laser system malfunction; recurrent corneal erosion; epithelial erosion at the incision site; epithelial defects; epithelial hyperplasia; superficial corneal opacification; incision opacification; endothelial perforation; incorrect entry of the channel; segment displacement; corneal melting; corneal edema; suspected infectious keratitis; infection following implantation (keratitis); microbial / bacterial keratitis; methicillin-resistant Staphylococcus aureus (MRSA) acute keratitis; late or early infection; sterile keratitis; inflammation around ICRS; glare; persistent photophobia; persistent foreign-body sensation; localized Descemet detachment; corneal thinning over the incision site; broken body of ring or segment broke; broken end orifice of ring; distant visual acuity worse than before surgery; ring tip near the incision site; sterile white/yellow deposits; progressive vision impairment or decreased vision; vacuum or suction loss during laser tunneling; subconjunctival hemorrhage due the suction ring of the dissector; re-intervention; segment override; residual haze in the area in which the explanted rings had been placed, epithelial ingrow in the ICRS tunnel.
- In general, these complications occurred in initial cases, resulting from the surgeon's learning curve or manual implantation technique. The cases of ring extrusion may be related to long arc rings, treatment with Cross-Linking performed simultaneously with the ring implant and patients with a history of atopy and eye friction. The safety and effectiveness of Keraring implantation have been established when used as directed. However, visual rehabilitation may be unsatisfactory if this device is not used in accordance with guidelines herein listed, by an experienced surgeon and with precision for dissection of the intrastromal tunnel, especially in the manual technique, or if patient fails to receive proper post-operative care.
- The safety and effectiveness of the Keraring implant are established when used according to the instructions for use, however, the visual result may be compromised if the device is not implanted according to the recommended surgical technique, by an experienced surgeon and with precision to dissection of the intrastromal tunnel, especially in the manual technique, or if the warnings and contraindications are not observed, or even if the patient does not receive adequate postoperative follow-up.

- Some patients with large, dilated pupil diameters are predisposed to low light visual symptoms and should be appropriately advised.
- It is anticipated that spectacle correction or contact lens correction may still be necessary following implantation of this device, especially in moderate to advanced keratoconus patients.
- Under mesopic conditions, patients may experience some loss in contrast sensitivity at low spatial frequencies.
- The safety and effectiveness of alternative refractive procedures following implantation of the Keraring has not been established.
- Some patients may have decrease of visual acuity after Keraring implantation. The patient should be informed about the possibility of visual acuity decreasing after the surgical procedure.
- Patients with higher levels of CDVA (Corrected Distance Visual Acuity) during preoperative or no vision loss after post-LASIK ectasia should not be considered suitable candidates for Keraring implant due to the risk of impaired vision after implantation.
- The progression of keratoconus is not related to the performance of the product, it is related to the evolution of the disease itself. The patient should be informed about keratoconus progression.
- The Keraring implant is reversible and adjustable for the treatment of keratoconus for each patient. In addition, in the case of segment explant after implantation of it without adequate success, good results can be achieved with the implant of a new segment, with significant visual and refractive improvements. The patient should be informed about the reversibility of implantation of Keraring.
- If there is any damage or defect in the product's label or instructions for use or any suspicion of quality deviation in the product, it should not be deployed and should be referred to Mediphacos Customer Service (SAC) so that the product is evaluated for its quality.
- If there is a suspicion of complications or adverse events associated with the Keraring implantation or serious incident with Keraring, Mediphacos' SAC, representative authorized and the competent authority of the Member State in which the users and/or patients are established must be informe.

11. PRECAUTIONS FOR USE AND STORAGE

- Do not use after sterilization expiration date indicated on package. The use of the expired product does not guarantee the maintenance of the sterility of the product.
- Single use only. Do not re-sterilize. The reprocess is forbidden. The reuse and/or re-disinfection and/or Re-sterilization of the product may not be effective and may cause inflammatory reactions, infection and contamination of the patient, toxic and mutagenic effects, hypersensitivity and/or explanation.
- The product should not be exposed to any solution other than sterile Balanced Saline Solution (BSS).
- This device must only be implanted by an experienced corneal surgeon trained in intrastromal corneal ring implantation techniques and management of related complications. Otherwise the intended use of the product may not be achieved.
- Handling of this device requires the use of appropriate techniques and instruments to ensure the integrity of the device and the performance of the surgical technique. The use of instruments provided by Mediphacos is recommended. The use of surgical instruments by a professional not trained in the ICRS implantation technique can cause corneal damage.
- Any device damaged during handling should not be implanted. The device should be discarded in the hospital waste according to the procedure established in the institution in which the surgical procedure would be performed.

- The correct selection of device dimensions and position for a specific patient depends on preoperative manifest refraction and the presence and extent of other corneal conditions, such as keratoconus. Refer to clinical guidelines and updated nomograms available from Mediphacos. Please, contact, keraring@mediphacos.com
- Do not use if sterile pouch is found opened and/or damaged. Damaged packaging can compromise product sterility, leading to risk of contamination of the patient. The device should be discarded in the hospital waste according to the procedure established in the institution in which the surgical procedure would be performed.
- Sterile packaging may have a sterilization indicator in beige and green, in addition to brown, to indicate that the product is sterile.
- Check the dimensions of the implant, the orientation (clockwise or counterclockwise) of the asymmetric rings and validity of the sterilization on the label before opening the sterile packaging.
- Store at room temperature (18 to 25 °C) and in the humidity range of 30 to 70 HR%. Product stored out of recommended conditions may impact its stability, compromising its intended use.
- If for any reason the device is removed from sterile condition and it is not implanted, then, the device must not be used and appropriately discarded on hospital waste according to the procedure established in the institution in which the surgical procedure would be performed.
- Implant identification labels supplied must be affixed to patient and/or hospital and/or surgeon files to in order to allow tracing of each individual device to its end user (patient)), if necessary, as required by national and international regulations. The implantation of the device should be performed in a sterile room with all the precautions associated with intraocular surgical procedures.

12. PATIENT PREPARATION

- If the patient wears contact lenses, he/she should be advised to discontinue its use at least two weeks before their final preoperative examination to allow accurate refraction and corneal topography measurements.
- Patients should be instructed not to wear makeup, especially around the eyes, on the day of surgery.
- Patients should be instructed not to rub their eyes after implantation to avoid migration or extrusion of the device and to ensure proper incision healing.
- Patients should be instructed on the importance of using all postoperative medication according to prescription.
- Patients should be instructed to contact their physician immediately if they experience pain, irritation, foreign body sensation, or if they experience a sudden change in vision after in the postoperative period.
- The patient should be instructed about the post-surgical follow-up to be performed to evaluate the surgical results achieved and complications, the periodicity of such follow-up should be defined by the surgeon.

NOTE: Failure to follow the care described above may compromise the results of the intended use of the product.

13. PATIENTS INFORMATION

- Implant card supplied must be filled by the surgeon and it must be given to the patient.
- Patients should be instructed that Keraring summary of safety and clinical performance is available in EUDAMED database: <https://ec.europa.eu/>. Until the Eudamed database is active, the Summary of Safety and Clinical Performance of Keraring will be available on the Mediphacos website (<https://mediphacos.com/produtos/cornea/keraring-anel-intraestromal/>).
- The possible extractable substances from the Keraring Intrastromal Corneal Ring Segment identified through

laboratory tests were Zinc and Bumetrizole (which is the UV blocker of Keraring - IUPAC name: 2-(2'-hydroxy-3'-tert-butyl-5'-methylphenyl)-5-chlorobenzotriazole)). According to the toxicological limit calculated for each substance, the amounts found are at least 6 million times lower than the toxicological limit.

14. GUIDELINES FOR SURGICAL TECHNIQUE

Perform cleaning, disinfection and sterilization of surgical instruments, according to the instructions of the instrument manufacturer. Before use, check that each surgical instrument is intact and without signs of degradation or damage, also check that they are completely clean, without any deposit, dirt or particle and sterile. Do not use surgical instruments if there are signs of degradation or damage or if the sterile package is open and/or damaged. Damaged packaging can compromise the sterility of the product, leading to the risk of contamination of the patient.

Open the sealed box and remove the envelope containing the sterile device. Open the sterile envelope and transfer the plastic case to sterile field. Remove the seal from the case. Place a few drops of sterile balanced salt solution (BSS) on the case lid and on the tip of the sterile tweezers, to eliminate any electrostatic charge and easy gripping of the implant by the tweezers. The case must be opened continuously and carefully, avoiding abrupt movements, to ensure the ring integrity. Using appropriate aseptic techniques, remove the segment from the kit using a specific forceps for manipulating intrastromal corneal rings. The use of the product without observation of the aseptic techniques can promote inflammatory reactions and contamination of the patient. Carefully inspect the device for any breakage and/or imperfections and/or adhered particles on the surface of the product to avoid inflammatory reactions and contamination of the patient. Proceed with the implant in the corneal stroma by specific surgical technique.

On the Mediphacos website can be found the surgical techniques recommended for the implantation of Keraring. The ring implant surgery takes a few minutes, varying according to the surgical technique, skill and experience of the surgeon.

MANUAL TECHNIQUE

In the mechanical or manual technique, the surgeon must mark the center of the pupil in order to use it as a reference point during the procedure. The pupillary center of the cornea is marked with a Sinskey hook. This is the reference point. For the manual procedure, a 1 mm long radial incision created 70-80% deep in the thickness of the curved meridian of the cornea, in the diameter of 5.0 mm or 6.00 mm, with the calibrated diamond scalpel. This is facilitated by intraoperative ultrasound pachymetry. Snap hooks are then used to create corneal pockets on either side of the radial incision. Then, two semicircular dis-sectors are introduced into the lamellar pouch by rotating clockwise and counterclockwise. Two semicircular tunnels are then created and the ring segments are inserted into the channels.

FEMTOSECOND LASER-ASSISTED TECHNIQUE

The femtosecond laser acts through the photodisruption mechanism promoting tissue dissection, forming micro bubbles of carbon dioxide and water vapor at the chosen depth. It has the benefit to create a tunnel with uniform and precise depth of 360 degrees, being therefore, safer and liable to less intraoperative complications than the manual technique.

Identification of the corneal optic center is done using a pen based on the corneal reflex of the surgical microscope coaxial light, asking the patient to fix the center of its filament. Then the patient is positioned under the laser and the suction and applanation system of the femtosecond laser is applied. The suction and applanation system are centered flattening the cornea and the laser is activated forming a tunnel of 360 degrees according to the chosen depth (usually 70% of the thickness of the thinner point in the 5-6 mm area, ring implant place)and with the incision on the corneal axis defined by the chosen Keraring type. Then, patency of the tunnel is ensured with a spatula passed through it followed by implantation of the Keraring segments.

The Keraring Intrastromal Corneal Ring Segment has an estimated lifetime of, at least, 3 years after implantation. After this period, the ophthalmologist should examine the patient, to verify that the ring continues to be safe and effective and that the ring's performance continues to meet its intended use.

15. WARRANTY AND LIMITATION OF LIABILITY

Mediphacos Indústria Médicas S/A guarantees that the product when delivered, is free from defect materials and in compliance with product's current version specifications. This warranty is valid until the expiration date of the product indicated on external labeling. Mediphacos may, in accordance with its policy, restore or arrange the exchange if any defects are found, as soon as the defective product returns to Mediphacos.

Mediphacos shall not be liable for any incidents and consequential loss and/or damage resulting directly or indirectly from the use or inability to use the product in any patient. Mediphacos neither assume, nor authorizes any person to assume on its name, any other liability responsibility with respect to this product other than set forth herein.

GAMINIO INFORMACIJA IR NAUDOJIMO INSTRUKCIJA - LT - REV15
Keraring Intrastromal Corneal Ring Segment | Modeliai: SIS, SI6, ASS, AS6, SG

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MDSS – Medical Device Safety Service GmbH